




GENCO
Tibbi Cihazlar



www.gencotibbi.com

SPONJEL



HAEMOSTATIC ABSORBABLE GELATINE SPONGE

- Rapidly stops the capillary and venous bleeding.
- Haemostatics time is 5 to 10 minutes.
- Ensures even coagulum.
- Provides clear and open surgical fields.
- Capable to absorb 50 times of its own weight.
- Absorbed by body within 3-4 weeks without producing antigenic reaction with thrombin.



SPONJEL GENERAL SIZES

80 mm x 50 mm x 10 mm

70 mm x 50 mm x 10 mm

80 mm x 30 mm (Anal)

10 mm x 10 mm x 10 mm (Dental)

Please contact us for other sizes

UNOWAX



BONE WAX

- UNOWAX is a haemostatic bone wax which is manufactured from natural beeswax.
- Enables to control the bleeding on the bone tissues via forming a mechanical barrier.
- Packed with double packaging which ensures 5 years of shelf life.



UNOWAX GENERAL SIZES

2.5 grams


3 grams

Please contact us for other sizes



PVA NASAL & EAR DRESSING

- Nasal, epistaxis, sinus, otology and other dressing sizes available.
- UNOPORE is made of PVA (Polyvinylalcohol) expansible sponge.
- Soft, comfortable to use.
- Vacuum sterile packaging ensuring long and safe protection.
- Micropore size minimizes tissue ingrowth.
- Capable to absorb 25 times of its own weight.
- Expansion capability under mucosal pressure.



UNOPORE GENERAL SIZES

80 mm x 15 mm x 20 mm


80 mm x 15 mm x 20 mm
(With Airway)

80 mm x 15 mm x 30 mm
(Anatomic Shape)

45 mm x 15 mm x 20 mm


45 mm x 15 mm x 20 mm
(With Airway)

Please contact us for other sizes



POLYPROPYLENE MESH

- UNOMESH is made of 100% nonabsorbable monofilament PP.
- Thickness of PP mesh is between 0.3 mm to 0.6 mm.
- Porosity of the mesh is 1.2 mm to 2.5 mm.
- UNOMESH knited with 80 μ to 160 μ filament to keep high quality.



UNOMESH GENERAL SIZES

6 cm x 11 cm

7.5 cm x 7.5 cm

10 cm x 15 cm

15 cm x 15 cm

30 cm x 30 cm

Please contact us for other sizes





UNOCEL & UNOCEL CA-L OXIDIZED CELLULOSE

- UNOCEL fully haemostats within 2-3 minutes.
- With UNOCEL's special high density texture, it easily controls large area bleeding surgery areas.
- UNOCEL absorbable in 12-72 hours and eliminated in 7-21 days.
- UNOCEL is bactericidal and bacteriostatic, kills wide range of gram positive and gram negative bacteria.

UNOCEL GENERAL SIZES

50 mm x 75 mm

50 mm x 350 mm

90 mm x 150 mm

100 mm x 200 mm

UNOCEL CA-L POWDER

2 grams

3 grams

4 grams

Please contact us for other sizes

* UNOCEL CA-L POWDER; manufactured from %100 oxidized cellulose with the active substance based on hydrogencalcium salt.



GENCO

Tıbbi Cihazlar

Genco Tıbbi Cihazlar San. Tic. Ltd. Şti.

1571 Sok. No: 22 Kat: 2 Çınarlı, İzmir / TÜRKİYE

Tel: +90 232 462 16 49 • Faks: +90 232 435 31 16

E-posta: info@gencotibbi.com