



Salutem —OZON—



Our Story;

As Has Medikal, we laid our foundations in Izmir about a quarter of a century ago. Since then, we continue to produce medical devices with determination that "help improve the health of people" by sticking to codes of conduct.

Our sole motivation is human health; and we are well aware of the responsibility it brings. This awareness pushes us to not only be a commercial brand, but also an effective supporter of the health industry. We are putting our efforts into becoming an international brand and to sustain our success.

Although our headquarters is located in Izmir, we have many branches next to State and Private hospitals and clinics from Edirne to Hakkari. We serve domestic and international markets with our strong partner "Salutem Ozon".

As we mentioned above, our motto is "health first", and we provide medical devices developed by experts according to this principle.

Stay well...













01

SALUTEM PLUS MEDICAL OZONE MACHINE

It is designed for professional use, institutional and university use, and all medical / dental / laboratory research applications. If you use ozone-resistant consumables while performing ozone application, you can increase the efficiency of the application.

Use methods;

Major autohemotherapy, minor autohemotherapy, rectal insufflation, topical therapy, saline (SF) etc. It is a full-featured ozone generator that allows you to implement all ozone protocols. It produces ozone concentration from 0.1 μ g/ml to 80 μ g/ml.

Technical Specifications

SCREEN	7 inch
LOW DOSE RANGE	0,1-0,2-0,3-0,4-0,5-0,6-0,7-0,8-0,9-1-1,1-1,2-1,3-1,4-1,5-1,6- 1,7-1,8-1,9-2-2,1-2,2-2,3-2,4-2,5-2,6-2,7-2,8-2,9
HIGH DOSE RANGE	3-4-5-6-7-8-9-10-11-12-13-14-15-16-17-18-19-20-21-22-23- 24-25-26-2771-72-73-74-75-76-77-78-79-80
PRESSURE SENSOR	YES
OXYGEN FLOW SENSOR	YES
VACUUM SENSOR	YES
NO OXYGEN WARNING	YES
NO OZONE WARNING	YES
VACUUM VALUE	150-200-250-300-350-400-500-600-700 mbar ayarlanabilir. 800 mbar sabit
TIMER	YES
CUPPING THERAPY VACUUM	150-200-250 mbar
COLOUR	Black
WARRANTY	1 years

^{*}Easily set up Sangiset, i-set settings.

^{*}Easily use for multi-plus applications.

^{*}Use the timer for physiological serum application, it automatically stops use when the timer ends and you will hear the sound alarm.

^{*}Stop ozone automatically with the timer during topical and rectal applications.







Vacuum citrate glass bottle



Transfusion seti



Transfer filtre kit



Butterfly needle



Injector



Rectal pump set



Topical application



Yenişehir Mah. 1145/1 Sok. No: 74 Konak / İzmir - TÜRKİYE

Tel: +90 232 433 38 33 Faks: +90 232 433 38 08

www.hasmedical.com / info@hasmedical.com



Parseller Mah. Taşköprü Cad. No: 33 - A Ümraniye / İstanbul Tel: +90 216 660 00 95 www.salutemozon.com info@salutemozon.com