

MaxFill® Bulk-Flow

BULK FLOWABLE COMPOSITE RESIN



DESCRIPTION

With the development of the bulk fill system, dentists have become to apply single-layer composite comfortably. This system provides not only low shrinkage stress but also effective light reflections for superior color harmony and aesthetics.

Specially designed glass filler materials provide naturalness with a chameleon effect on dentin and enamel with natural light transmission and reflection.

INDICATIONS

- Maxfill Bulk Flowable can be easily and precisely applied from the syringe. For class I and II cavities it is suitable as a baseliner (pedestal material) or as a lining and filling material for small posterior cavities.
- It is an excellent carcass production material and has high occlusal strength. Its smooth, non-sticky structure also makes it ideal for direct posterior restorations.
- Single-stage application, does not require an additional top-end layer.
- It offers excellent adaptability, manipulation, and formability without the need for expensive application devices.

ADVANTAGES

- Low polymerization shrinkage and shrinkage stress
- Excellent hardening depth up to 4 mm thick
- The most accurate color matching (Chameleon effect)
- High radiopacity



Discover
the potential

