



# RESTORING SMILES ONE PRODUCT AT A TIME

PRODUCT CATALOGUE  
2022 - 2023

[www.saridental.com.tr](http://www.saridental.com.tr)



# INDEX

## INTRODUCTION

### RESTORATIVES

Crown Pit and Fissure Sealant	5
Crown Glass Liner	5
CrownFill Nano Composite	6
CrownFill P Posterior Composite	6
CrownFlow Composite	7
CrownEtch Gel	7
CrownBond	8
CrownSE Bond	8

### ENDODONTICS

CrownEdta Gel	11
Crown Pex	11
CrownPaste	12
CrownPaste LC	12
CrownFile Rotary Files	14
CrownFile Reciproc Files	14
CrownFile Pedo Files	15
CrownFile Path File	15

### PROPHYLAXIS

Crown Gingiva Barrier	21
-----------------------	----

### POLISHING & FINISHING

Crown Polishing Discs	23
-----------------------	----

### AUXILIARY

CrownRemover	25
--------------	----



# ABOUT US

**SARI DENTAL** WAS ESTABLISHED IN 2005 IN ISTANBUL. **SARI DENTAL** HAS BEEN SELLING, PRODUCING, EXPORTING AND IMPORTING DENTAL CONSUMABLES SINCE ITS ESTABLISHMENT.

THE PRODUCTION OF MANY MATERIALS SUCH AS DENTAL AMALGAM FILLING MATERIAL AND DENTAL COMPOSITE FILLING MATERIAL UNDER THE CROWNAL-LOY BRAND IS CARRIED OUT BY OUR COMPANY. IN ADDITION TO OUR DOMESTIC SALES, IT EXPORTS DENTAL CONSUMABLES TO MORE THAN 5 COUNTRIES.

BEHIND THE SUCCESS OF **SARI DENTAL** IS A CUSTOMER-ORIENTED AND PRODUCTIVITY-CENTERED MANAGEMENT APPROACH. THIS APPROACH DEVELOPS NOT ONLY WITH MATERIAL GAINS, BUT ALSO WITH THE AWARENESS OF CORPORATE CITIZENSHIP THAT THE WHOLE SOCIETY BENEFITS FROM AND WILL BENEFIT FROM.

WE HAVE REGISTERED OUR UNDERSTANDING OF QUALITY WITH ISO 13485 AND CE CERTIFICATES.

IN THIS CONTEXT, IT IS AMONG OUR AIM TO ENSURE THAT ALL THE PRODUCTS WE SELL ARE THE FIRST BRANDS THAT COME TO MIND FOR OUR DENTISTS.

## **SARI DENTAL, IN THE DENTAL INDUSTRY;**

MAKING THE DIFFERENCE,

RESEARCHING, FINDING AND USING THE RIGHT TECHNOLOGIES AT HOME AND ABROAD

QUICKLY ADAPTING TO TECHNOLOGICAL DEVELOPMENTS

DRIVING CHANGE

IT IS A COMPANY THAT MAKES QUALITY PRODUCTION AND KEEPS USER AND EMPLOYEE SATISFACTION AT THE HIGHEST LEVEL. IT AIMS TO BE AN EXAMPLE IN EVERY ASPECT OF ITS PRODUCT, SERVICE AND TRADE.







# Crown Pit and Fissure Sealant

Light curing sealant with fluoride releasing for occlusal surfaces. It can be used as preventive products.

## PROPERTIES

Fluoride releasing  
Radiopaque  
Easy application  
Low shrinkage

## PACKAGING

2 gr syringe  
2 gr syringe x 4



# Crown Glass Liner

Glass Liner is a radiopaque and light-cured cavity sealant with high strength and good adhesion to dentin. Glass Liner contains glass ionomer and fluoride, it is applied directly. It is used under all filled materials.

## PROPERTIES

Durability: 200 MPa  
Hardening Depth: 5mm

## PACKAGING

2 gr syringe  
2 gr syringe x2  
2 gr syringe x4  
Syringe Tips



# CrownFILL

CrownFill Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

## PROPERTIES

Universal Nano Composite It has light refraction property like a normal tooth due to the nano fillers in it It has high filler rate It can be polished high with high degree It has high resistance to abrasion It provides ease of use

## PACKAGING

Refill – 4 gr in a syringe

Kit 1 - 4 gr syringe x 6 pcs. , 5 ml CrownBond x 1 pcs. 3 ml CrownEtch x 1 pcs. and Accessories

Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach



# CrownFill Posterior

CrownFill Posterior is micro hybrid composite system which can be the photocurable applied in the posterior field.

## PROPERTIES

Non stick to instruments It has surface hardness is like real enamel It has high intensity It has a perfect consistence can be shaped/formed easily It has high filler rate It resists to chewing force There is not discoloration after the polymerization.

## PACKAGING

Refill – 4 gr in a syringe

Kit 1 - 4 gr syringe x 6 pcs. , 5 ml Crown Bond x 1 pcs. 3ml CrownTech x 1 pcs. and Accessories

Shades – A1 – A2 – A3 – A3,5



# CrownFlow

CrownFlow is microhybrid flowable light curing radiopaque composite.

## PROPERTIES

Low shrinkage, low wear High gloss polishability Beautiful esthetic restorations Very low surface roughness Excellent handling Non sticky Very good marginal adaptation

## PACKAGING

2 gr in a syringe  
2 gr syringe x 4 pcs.  
2 gr syringe x 4 pcs. , 5 ml Crown Bond x 1 pcs. 3ml Crown-Tech x 1 pcs. and Accessories Shades – A1 – A2 – A3 – A3,5 – B1 – OA2



# CrownEtch

CrownEtch is an etching gel containing %37 phosphoric acid. Delivered in 3 ml. and 50 ml syringe with very thin application needles.

## PROPERTIES

Easily Rinse  
Blue coloured  
Perfect Surface Retention

## PACKAGING

3 ml syringe  
3 ml syringe x 4  
50 ml syringe

## CrownBond

CrownBond is light curing, single component bonding agent for the adhesive restorations. It universal and compatible with all composites.

### PROPERTIES

Low abrasion Elastic connection

### PACKAGING

5 ml in a bottle 5 ml in a bottle



## CrownBond SE

CrownSEBond is self-etching, light curing one component adhesive for the adhesive restorations.

### PROPERTIES

Very strong bonding to enamel and dentine No extra etching with phosphoric acid Long lasting restorations Time savings


### PACKAGING

5 ml in a bottle









# Endodontics



## Crown Pex

Preparation gel for root canal cleaning and lubrication

### PROPERTIES

Easy mobility  
Antibacterial effect  
Reduces file breakage

### PACKAGING

2gr syringe



## Crown Edta Gel

It is root canal preparation gel for lubricating and cleaning of root canal. It comes with application tips.

### PROPERTIES

Easily moving  
Antibacterial effect  
Reduce file breakage

### PACKAGING

3 ml syringe



## Crown Paste

Calcium hydroxide premixed paste non-setting for root canal therapy and vital pulp therapy. It comes with application tips.

### PROPERTIES

High Radiopacity  
Easy removal  
Paste has pH > 12  
Excellent handling characteristics

### PACKAGING

2 gr syringe



## Crown Paste LC

Light curing calcium hydroxide for pulp capping and liner material under dental filling.

### PROPERTIES


Direct & indirect pulp capping  
Light curing  
Single component

### PACKAGING

2 gr syringe











# CROWNFILE











HEAT TREATED CONTROL  
MEMORY ROTARY FILES



SX	S1	S2	F1	F2	F3	F4	F5
							
Starter	Shaping		Finishing				
Torque: 3.0 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

Excellent Cutting Efficiency

SX	S1	S2	F1	F2	F3	F4	F5
							
Starter	Shaping		Finishing				
Torque: 3.0 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

Maximum Flexibility

CROWNFILE Rotary Files manufactured using highest quality nickel titanium alloy, designed special and innovative heat treatment production process which gain properties of Control Memory with excellent cutting efficiency and great "exibility" to the instruments. After the special heat treatment production process, the wire becomes Gold or Blue colours. The innovative heat treatment production process has upgraded metallurgy exclusively developed by Sari dental.

Action Kit SX - S1 - S2 - F1 - F2 - F3

### Action Gold Rotary Files

Rotation Files  
Variable Taper  
Better Cutting Efficiency  
x2 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank

Wire : Gold  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section



### Action Blue Rotary Files

Rotation Files  
Variable Taper  
Better Flexibility  
X4 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank

Wire : Blue  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section



	10/06	15/06	20/06	25/06	30/06	35/06	40/06	45/06	50/06
	Shaping		Finishing						
Torque: 3.0 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

Excellent Cutting Efficiency

Maximum Flexibility

	10/06	15/06	20/06	25/06	30/06	35/06	40/06	45/06	50/06
	Shaping		Finishing						
Torque: 3.0 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm

CROWNFILE Rotary Files manufactured using highest quality nickel titanium alloy, designed special and innovative heat treatment production process which gain properties of Control Memory with excellent cutting efficiency and great "exibility" to the instruments. After the special heat treatment production process, the wire becomes Gold or Blue colours. The innovative heat treatment production process has upgraded metallurgy exclusively developed by Sari dental.

Touch Kit 15/06-20/06-25/06-30/06-35/06-40/06 19/08-15/04-20/04-20/06-25/06-30/06 19/08-10/06-15/06-20/06-25/06-30/06 19/08-15/06-20/06-25/06-30/06-35/06

### Touch Gold Rotary Files

Rotation Files  
Regular Taper  
Better Cutting Efficiency  
x2 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank

Wire : Gold  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section



### Touch Blue Rotary Files



















Rotation Files  
Regular Taper  
Better Flexibility  
X4 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank

Wire : Blue  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section






	10/04	15/04	20/04	25/04	30/04	35/04	40/04	45/04	50/04
									
									
<b>04 TAPER</b>	Starter		Shaping / Finishing		Finishing				
	Torque: 3.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm

Excellent Cutting Efficiency

Maximum Flexibility

	10/04	15/04	20/04	25/04	30/04	35/04	40/04	45/04	50/04
									
									
<b>04 TAPER</b>	Starter		Shaping / Finishing		Finishing				
	Torque: 3.0 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm	Torque: 1.5 Ncm 300 Rpm

CROWNFILE Rotary Files manufactured using highest quality nickel titanium alloy, designed special and innovative heat treatment production process which gain properties of Control Memory with excellent cutting efficiency and great "exibility" to the instruments. After the special heat treatment production process, the wire becomes Gold or Blue colours. The innovative heat treatment production process has upgraded metallurgy exclusively developed by Sari dental.

SMART Kit 15/04-20/04-25/04-30/04-35/04-40/04 19/08-10/04-15/04-20/04-25/04-30/04 19/08-15/04-20/04-25/04-30/04-35/04 19/08-16/02-20/04-25/04-30/04-35/04

## Smart Gold Rotary Files

Rotation Files  
Regular Taper  
Better Cutting Efficiency  
x2 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank

Wire : Gold  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section





## Smart Blue Rotary Files

Rotation Files  
Regular Taper  
Better Flexibility  
X4 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank

Wire : Blue  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section





CROWNFILE Rotary Files manufactured using highest quality nickel titanium alloy, designed special and innovative heat treatment production process which gain properties of Control Memory with excellent cutting efficiency and great flexibility to the instruments. After the special heat treatment production process, the wire becomes Gold or Blue colours. The innovative heat treatment production process has upgraded metallurgy exclusively developed by Inci Dental.

Pedo Kit 15/06-20/04-25/04 15/06-20/04-30/04

## PEDO Gold Rotary Files

Rotation Files  
Specially produced for pedodontic cases  
Regular Taper  
Better Cutting Efficiency  
x2 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Taper and Size information on the shank

Wire : Gold  
Length : 18 mm  
Packaging : 3 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section



## PEDO Blue Rotary Files


Rotation Files  
Specially produced for pedodontic cases  
Regular Taper  
Better Flexibility  
X4 times better Cyclic Fatigue Resistance than conventional Files  
Excellent Strength  
No bounce back  
Respect Original Anatomy of Root Canal  
Great Fracture Resistance  
Lasered Taper and Size information on the shank

Wire : Blue  
Length : 18 mm  
Packaging : 3 Files

Control Memory  
Heat Treated Wire  
Convex Triangle  
Cross - Section







Control Memory  
Heat Treated Wire

Convex Triangle  
Cross - Section




Wire : Gold  
Length : 18 mm  
Packaging : 3 Files

Retreatment Kit D1-D1-D2-D2-D3-D3

### Rotation Files

Specially produced for retreatment cases  
Designed to remove filling materials like Gutta Percha  
x2 times better Cyclic Fatigue Resistance than conventional Files  
Better Cutting Efficiency  
Great Fracture Resistance  
Lasered Taper and Size information on the shank



Control Memory  
Heat Treated Wire

Convex Triangle  
Cross - Section



Wire : Blue  
Length : 18 mm  
Packaging : 3 Files

### Rotation Files

Specially produced for creating glide path  
Easily reach apex  
Better Flexibility  
x2 times better Cyclic Fatigue Resistance than conventional Files  
Great Fracture Resistance  
Lasered Length Marked  
Lasered Taper and Size information on the shank







# PROPHYLAXIS





# PROPHYLAXIS

## Crown Gingiva Barrier

Light curing resin-based gingiva barrier to protect soft tissue.

### PROPERTIES


- Easily removable after treatment
- It minimizes heat and soft tissue sensitivity during curing
- Blue shade
- 5 second curing

### PACKAGING

2 gr syringe







# POLISHING & FINISHING

# CrownDisc Polishing Discs

Has different kind of discs which has different grins and different sizes.  
It allows excellent results in polishing.  
Flexible and easy to use.

## PROPERTIES

10 mm blue, yellow, white, green, brown

12,7 mm blue, yellow, white, green, brown

14 mm blue, yellow, white, green, brown

White – Extra Fine for Polishing

Yellow – Fine for Finishing

Green – Medium for Countoring

Blue – Coarse for Gross Reduction

Brown – Extra Coarse for Gross Reduction


## PACKAGING

40 pcs discs and 1 mandren in a plastic box

100 pcs discs and 1 mandren in a plastic box







# AUXILIARY

# Crown Remover


It is practical crown remover for easy procedure. It removes crown using the way of impact method with arc.

## PROPERTIES

Can be sterilized at 134 °C

3 types of remover tips

Different types of impact





# SARI DENTAL

[www.saridental.com.tr](http://www.saridental.com.tr)

E-mail: [info@saridental.com.tr](mailto:info@saridental.com.tr)

Phone : 0212 523 58 71 - +90 542 309 55 40

Address: Molla Gürani Mah. Oğuzhan Cad.  
Oğuzhan Apt. No:19/D Fındıkzade - Fatih / İSTANBUL



AVAILABLE ON



WEBSITE



@saridental.official