

**KAFGRUP**

**ALCOHOL-BASED HAND AND SKIN ANTISEPTIC/DISINFECTANT GEL/  
DISINFECTANT WIPES AND SWAB WIPES**

**ALCOHOL-BASED/FREE RAPID SURFACE DISINFECTANT WIPES AND  
CONCENTRATED FLOOR/ONESURFACE WIPES**

**ANTISEPTIC SOLUTION/ANTIBACTERIAL MEDICAL SOAP**

**ULTRASOUND ECG GEL/IPL GEL**

**ANTIBACTERIAL,ANTIFUNGUSIDAL AND ANTIVIRAL BODY CLEANSING WIPES/  
INCONTINENCE BODY CARE PRODUCTS**

**WET WIPES / WET BABY TOWEL**

**SODIUM HYPOCHLORIC/EDTA SOLUTION**

**ADHESIVE REMOVER SPRAY / EEG PASTE**

## ALCOHOL-BASED HAND AND SKIN ANTISEPTIC/ DISINFECTANT GEL



Gel/Liquid	50 ML
Gel/Liquid	100 ML
Gel/Liquid	250 ML
Gel/Liquid	500 ML
Gel/Liquid	1 L
Gel/Liquid	4 L
Gel/Liquid	5 L

- OneHand is an alcohol-based hand and skin antiseptic.
- Used for hygienic and surgical hand disinfection.
- Does not cause any allergic reaction.



Antivirus



Kills 99.9%  
Of tested  
bacteria  
fungi and  
viruses



Ethanol

## DISINFECTANT WIPES AND SWAB WIPES



Disinfectant Wipes	1 PCS
Alcohol Swab	1 PCS
Chlorhexidine Swab Wipes	1 PCS

- Ideal for fast and effective cleaning, WanCare Disinfectant Wipes very useful and practical for personal hygiene and kills 99.9% of bacteria and viruses.



Antivirus



Kills 99.9%  
Of tested  
bacteria  
fungi and  
viruses



Ethanol

## ALCOHOL-BASED/FREE RAPID SURFACE DISINFECTANT WIPES AND CONCENTRATED FLOOR - SURFACE DISINFECTANT



Spray Bottle	750 ML
Spray Bottle	1 L
Wipes	100 PCS
Wipes	200 PCS
Wipes	350 PCS
Canister	5 L



Antivirus



Kills 99.9%  
Of tested  
bacteria  
fungi and  
viruses

- It is a fast surface disinfectant used on floors and surfaces.

## ONESURFACE WIPES



KAF G36-1	100 PCS
KAF G36-2	150 PCS
KAF G36-3	200 PCS

- Usage Area:** OneSurface Wipes is used for rapid disinfection of non-alcoholic and non-invasive medical devices used for rapid disinfection and cleaning on medical and medical device surfaces, incubators, aerators used in dentistry and medical devices with obligatory disinfection used in ENT examinations and materials that are not resistant to alcohol.



Alcohol,  
aldehyde and  
phenol free



Ready to use



Fast and effective  
disinfection



Sporicidal  
activity in 1  
minute



## ANTISEPTIC SOLUTION



PI 4

1 L - 10%

PI 15

1 L - 7,5%

7,5%

Povidone  
Iodine

10%

Povidone  
Iodine

- **Usage Area:** It is used in hygienic and surgical hand disinfection, in the disinfection and preparation of the surgical site before the operation, in the disinfection of the skin and around the wound, in antiseptic washings.

## ANTIBACTERIAL MEDICAL SOAP



KAF G72

1 L

KAF G72-1

5 L

4%

Chlorhexidine

- **Usage Area:** Antibacterial medical soap is an antiseptic, antimicrobial skin cleanser used by medical professionals before surgical procedures and by patients before a surgical procedure. This special soap cleans the surgeon's own skin as well as their patients. It is essential that all germs are eradicated on the skin surfaces of surgeons, operating room attendants, and the patient in order to help prevent postoperative infections.

## ULTRASOUND ECG GEL



Blue		Transparent
KAF G31-4	250 ML	KAF G31-2
KAF G31-5	500 ML	KAF G31-1
KAF G31-6	1 L	KAF G31
KAF G31-7	5 L	KAF G31-3



Hypoallergenic



Odorless



Easy to Clean



100%

Odorless

## Ultrasound Ecg Gel/Ipl Gel

- **Usage Area:** Ensures clearer and uninterrupted transmission of ultrasonographic waves by disconnecting the air between the skin and the probe during any and all ultrasonography, doppler ECG or effort test. The active substances contained do not damage the probe thanks to the high viscosity and the vacuum technology employed for production thereof. During operation, does not become affected by the body salt, ensures comfortable ultrasonography. Recommended for use by a specialist MD.
- **Ingredients:** Deionized water, Carbomer, TriEthanol Amine, MonoPropylene Glycol, Benzyl Alcohol, Methylchloroisothiazolinone

## IPL GEL



Transparent	
K03	250 ML
K03-01	500 ML
K03-02	1LT
K03-03	5LT



Hypoallergenic



Odorless



Easy to Clean



100%

Odorless

## ANTIBACTERIAL, ANTIFUNGUSIDAL AND ANTIVIRAL BODY CLEANSING WIPES



KAF G42

50 PCS

KAF G42-1

40 PCS

Kills Covid-19  
Virus

1 min



pH 5.5

- Total body cleansing with anti-bacterial effect
- No need extra dry and rinse after use
- Nourishes and moisturizes the skin with vitamin E, vitamin B, aloe vera and phospholipids.
- Protects against microorganisms with Sodium Benzoate and Potassium Sorbate Cleans, Soften and Soothes
- Prevents bacteria with Chlorhexidine Digluconate

## INCONTINENCE BODY CARE PRODUCTS



3%

Contains 3%  
DimethiconeBarrier  
effectCleans, softens  
and soothesLatex, lanolin,  
alcohol, Paraben  
freeNo rinsing  
or drying  
required

- **Usage Area:** Incontinence Body Care Products makes a barrier protective layer in order to protect the skin of patients with special needs, such as those with limited mobility, bedridden and intensive care patients, from external factors. it provides complete body cleaning and care for patients who cannot bath.
- **Composition:** Aqua, Phospholipid, Dimethicone, Glycerin, EDTA, Citric Acid, Polysorbate 20, Phenoxyethanol, Caprylyl Glucol.



## WET WIPES



- Provides effective cleaning on the skin thanks to its special formula. Formulated with skin compatible pH range. Dermatologically tested. Cleans, soothes and moisturizes the skin gently.



KAF G48-1	40 PCS
KAF G48	50 PCS
KAF G48-2	60 PCS
KAF G48-3	72 PCS
KAF G48-4	90 PCS
KAF G48-5	100 PCS
KAF G48-6	120 PCS



Paraben



Daily Care



Balanced

Dermatologically  
Tested

## WET BABY TOWEL



- Provides effective cleaning on the skin thanks to its special formula. Formulated with skin compatible pH range.



KAF G62-1	40 PCS
KAF G62	50 PCS
KAF G62-2	60 PCS
KAF G62-3	72 PCS
KAF G62-4	90 PCS
KAF G62-5	100 PCS
KAF G62-6	120 PCS



Paraben



Daily Care



Balanced

Dermatologically  
Tested



## SODIUM HYPOCHLORIC



D1

250 ML

5%

Sodium  
Hypochloric

- **Intended Use:** Eliminates soft and hard tissue residues from endodontic instruments and equipment.
- **Contents:** Sodium Hypochlorite, Distilled Water

## EDTA SOLUTION



D2

250 ml - 5%

D2-1

250 ml - 7%

D2-2

250 ml - 15%

5%

Edta  
Solution

7%

Edta  
Solution

15%


Edta  
Solution

- **Intended Use:** Cleans dental surfaces. In endodontic treatment, it is mainly used to facilitate canal instrumentation in attrition canals.
- **Contents:** Ethylene Diamine Tetra Acetic Acid, Sodium Salt, Sodium Hydroxide, Distilled Water

## ADHESIVE REMOVER SPRAY



Spray	ONCR002	50 ml
Spray	D4	100 ml
Aerosol	D5-1	100 ml
Wipes	D6	1 pcs



- Provides atraumatic application.
- It is effective in instantly removing the strongest adhesive and adhesive residues.
- It eliminates the side effects of medical adhesives on the skin.
- It is used to facilitate the painless and trauma-free removal of adhesive dressings, plasters and other adhesive products from the skin in a few seconds. It can be used safely on the skin of children, adults and the elderly. The remover is squeezed all over the product, waited for a few seconds and the adhesive is gently pulled from the skin.

## EEG PASTE



KAF G55	400 GR
KAF G56	250 GR

- The adhesive quality of OneCare EEG Paste is usually sufficient to adhere the electrodes to the skin and provide conductivity for the neurodiagnostic exam. Place the amount of paste needed on a surface such as a strip of surgical tape. Gently abrade the skin. Using a bell-shaped cup electrode with the hole in the center, place paste inside the electrode cup until it is slightly overfilled. Place the electrode onto the prepared electrode site and press with medium pressure. Use either tape or a piece of gauze to fix the electrode in place.



# YOUR PRIVATE LABEL!



## You may Introduce the Wancare Products with Your Own Brand ...

You may offer Wancare products in your own shops with your own brand. We also produce our products without a logo with the aim of bringing Wancare brand to large masses. You can create a good option for your customers and move your brand to better points with the high quality Wancare products.

KAF GRUP SAĞLIK HİZMETLERİ İNŞ. SAN. VE TİC. LTD. ŞTİ.

T +90 212 471 4 200 (pbx) / F + 90 212 471 4 201

Atakent Mah. 221. Sk. No:3A Rota Office A Blok Kat:14 D:83 K.ÇEKMECE / İSTANBUL

www.kafgrup.com www.wancare.com.tr info@kafgrup.com @wan.care & kafgrup

# KAF GRUP